


Impala 
ventilator

LifeKit by **MTTS**



IMPALA

Compact ventilator suitable for ambulances and intensive care wards




Overview

- Portable
- Suitable for emergencies or intensive care wards
- Build-in battery - 4 hours backup
- Integrated blender
- FiO₂ monitoring
- Easy-to-use touch screen
- ACV, SIMV and CPAP modes
- Controlled target PIP and consistent PEEP
- No time restriction on use due to the gas tank limitation

Reusable breathing circuits available



User friendly control panel



Physical Specifications (excluding stand)

Unit dimensions	20 cm x 22 cm x 20cm
Shipping (boxed) dimensions	30 cm x 25 cm x 30 cm
Net weight	3 kg
Shipping weight	4 kg

Performance Specifications

Power characteristics	100 - 240 VAC, 50-60 Hz
Integrated air compressor	Blower
Delivered oxygen concentration	21 - 100%
Battery	Li-ion, 12V-6000mAh

Ventilation mode	ACV, SIMV, CPAP
Respiratory rate	Adult: 5 - 40 bpm Neonatal/Pediatric: 5 - 120 bpm
Inspiration time	Adult: 0.1 - 6 s Neonatal/Pediatric: 0.4 - 6 s
Tidal volume	Adult: 5 - 800 ml Neonatal/Pediatric: 5 - 300 ml
PIP	2 - 40 cm H ₂ O
Pressure limitation	4 - 45 cm H ₂ O
PEEP	0 - 25 cm H ₂ O

Environmental Specifications

Operating	Ambient temperature: 0°C to +40°C Humidity: 0% to 90% RH non-condensing
Transport and Storage	Ambient temperature: -20°C to +60°C Humidity: 0% to 90% RH non-condensing

Standard Compliance

EN ISO 13485:2016	EN 60601-1-8:2007/A11:2017
EN ISO 14971:2012	ISO 80601-2-12:2020
EN 60601-1:2006/A1:2013	MEDDEV 2.12-1 Rev. 8
EN 60601-1-2:2015	MEDDEV. 2.7.1 Rev. 4
EN 60601-1-6:2010	MEDDEV. 2.12-2 Rev. 2

Ordering Information and Servicing

MAIN UNIT	ITEM ORDER NO
MTTS Impala Ventilator	VEN-10000
Reusable Circuit - adult	VEN-10030
Single use Circuit - adult	VEN-10040
Single use HME-HEPA filter	VEN-10050
Single use 3 in 1 filter	VEN-10051
Single use swivel straight connector	VEN-10052
Silicon mask - size 3	VEN-10063
Silicon mask - size 4	VEN-10064
Silicon mask - size 5	VEN-10065
Oxygen Hose	VEN-10070
Stand	VEN-10080
Test Lungs	VEN-10090

WARRANTY		
Standard Warranty	1 Year	Subject to MTTs policy

MTTS

MEDICAL TECHNOLOGY TRANSFER AND SERVICES Co., LTD
ISO 13485: 2016 Certified Company

For more information or to place an order please contact us at:

@ info@mtts-asia.com

www http://mtts-asia.com

+84 224 37 666 521

No 26, Lane 41 An Duong Vuong Street
Phu Thuong Ward, Tay Ho District
Hanoi, Vietnam

Distributor: