

Wallaby[™]
warmer

LifeKit by **MTTS**



CE 2265

WALLABY

An efficient, comfortable device to **prevent and treat hypothermia.**



Overview

- REMOTE TEMPERATURE and object sensing providing full patient protection from overheating if standard skin sensor is disconnected or faulty
- Angle adjusted bed (up to 12 degrees)
- Adjustable overhead unit and tray for X-ray procedures

1. Temperature Control System

- Automatic (servo controlled)
- Manual
- Prewarm

2. Displays/Indicators

- bright LCD display to monitor patient temperature, treatment time, total usage time, heater power level and sensor connection
- built in APGAR score timer
- LCD display with diagnostic and troubleshooting information
- total system time counter (lifetime)

3. Safety

- Heater element failure
- Sensor failure
- Power failure
- Sensor off the patient

4. Accessories

- IV Pole
- Storage Shelf
- Storage Drawer
- Cylinders holder



Dimensions and Weight

Unit dimensions	64 cm x 54 cm x 185 cm
Shipping (boxed) dimensions	38 cm x 82 cm x 122 cm
Net Weight	64 kg
Shipping Weight	76 kg

Performance Specifications

Application	Treatment of neonatal and pediatric patients to prevent hypothermia
Heater Element	Type: Radiant Power: 650W
Procedure Lamp	Type: LED Power: 3W

Temperature

Temperature range	34.0°C - 37.5°C in 0.1°C
Temperature measurement	10.0°C - 50.0°C in 0.1°C

Displays / Indicators

Set temperature	Total usage time
Patient temperature	Heater power level
Temperature Alarm	Sensor connected
Treatment time	Power indicator
APGAR score timer	

Extra Features

- Adjustable overhead unit for X-ray procedures
- Tray for X-ray films
- Adjustable bed angle (up to 12 degrees)
- Shelf for monitoring equipment

Requirements

Power Requirements	Voltage: 220V AC 10% Frequency: 47 - 63 Hz Wattage: 650W maximum
Agency Requirements	EN ISO 13485: 2012 EN 60601-1:2006/AC2010 EN 60601-1-2:2007/AC2010 EN 60601-6:2010 IEC 60601-2-21: 2003

Ordering Information and Servicing

MAIN UNIT	ITEM ORDER NO	
Warmer Unit	WLB-1000	Main Unit
ACCESSORIES	ITEM ORDER NO	
Skin Sensor - tape	WLB-1011	
Sticky tape for skin sensor (100)	WLB-1013	
Mattress	WLB-1020	
IV Pole	WLB-1030	
Storage Shelf	WLB-1040	
Storage Drawer	WLB-1050	
Cylinders Holder	WLB-1060	
Baby Bonnet	WLB-1070	
Documents Holder	WLB-1080	

WARRANTY

Standard Warranty	1 year	Subject to MTTs policy
-------------------	--------	------------------------

MTTS

MEDICAL TECHNOLOGY TRANSFER AND SERVICES Co., LTD
ISO 13485: 2016 Certified Company

@ info@mtts-asia.com

www http://mtts-asia.com

+84 243 766 6521

No 26, Lane 41 An Duong Vuong Street
Phu Thuong Ward, Tay Ho District
Hanoi, Vietnam

Distributor: