

Beluga   
resuscitator

LifeKit by **MTTS**



CE 2265


# BELUGA

An easy-to-use t-piece resuscitator suitable for emergencies and delivery units



## Overview

- Portable
- 4-hour battery backup
- Integrated blender
- Oxygen monitoring
- Easy-to-use touch screen
- Controlled target PIP and consistent PEEP
- No time restriction on use due to the gas tank limitation
- Usable in horizontal or vertical position



## Physical Specifications

Unit dimensions	20 cm x 22 cm x 20 cm
Shipping (boxed) dimensions	30 cm x 30 cm x 40 cm
Net weight	2.5 kg
Shipping weight	4.0 kg

## Performance Specifications

Power characteristics	100 - 240VAC, 47/63Hz
Manometer accuracy	± 1 cm H <sub>2</sub> O
Peak Inspiratory Pressure (PIP)	5 - 45 cm H <sub>2</sub> O
Positive End Expiratory Pressure (PEEP)	1 - 20 cm H <sub>2</sub> O
Delivered oxygen concentration	21 - 100%
Integrated air compressor	Blower
Battery	Li-ion, 12V-4000mAh

## Environmental Specifications

Operating	Ambient temperature: 18°C - 35°C Humidity: 0% to 90% RH non-condensing
Transport and Storage	Ambient temperature: 0°C - 50°C Humidity: 0% to 90% RH non-condensing

## Standard Compliance

EN 60601-1:2006/A1:2013  
EN 60601-1-2:2015

## Ordering Information and Servicing

MAIN UNIT	ITEM ORDER NO
MTTS Beluga Resuscitator	RES-10000
Reusable tube	RES-10011
Single use PEEP valve with T-piece	RES-10012
Single use T-piece circuit*	RES-10020
Silicon mask - size S (42mm)	RES-10060
Silicon mask - size M (50mm)	RES-10061
Silicon mask - size L (60mm)	RES-10062
Oxygen Hose	RES-10070
Arm Bag	RES-10080
Power cable - 12V Car Lighter*	COR-10006
*Optional	
WARRANTY	
Standard Warranty	1 Year Subject to MTTs policy

# MTTS

MEDICAL TECHNOLOGY TRANSFER AND SERVICES Co., LTD  
ISO 13485: 2016 Certified Company

For more information or to place an order please contact us at:

@ info@mtts-asia.com

www http://mtts-asia.com

+84 24 37 666 521

No 26, Lane 41 An Duong Vuong Street  
Phu Thuong Ward, Tay Ho District  
Hanoi, Vietnam

Distributor: