


Koala[®]
infant bed


LifeKit by **MTTS**



CE

KOALA INFANT BED

A safe, comfortable **bed** for **neonatal** and **paediatric** patients.



Overall

- suitable for neonates and paediatric patients
- compatible with MTTS Overhead Phototherapy
- rigid, steel body
- biocompatible materials
- safe and comfortable

1. Bassinet


- clear polycarbonate for easy patient observation
- removable side panels
- comfortable handhold around the rim
- smooth surface for quick wipe-down

2. Mattress

- waterproof fabric
- biocompatible material

3. Body - IV Pole

- adjustable height
- stainless steel



Dimensions and Weight

Unit dimensions	96.2 cm x 59 cm x 98.5 cm
Box dimensions	40 cm x 90 cm x 65 cm
Unit weight	44.0 kg
Shipping weight	55.0 kg
Mattress dimensions	54 cm x 64 cm x 3 cm
Outer bed dimensions	59 cm x 70 cm x 21 cm
Bed height	75 cm

Ordering Information and Servicing

MAIN UNIT	ITEM ORDER NO	
Infant Bed	IBD-4000	Main Unit
ACCESSORIES	ITEM ORDER NO	
Mattress	WLB-1020	
Wheel	IBD-4010	
Side Panel (L)	IBD-4021	
Side Panel (S)	IBD-4022	
X-ray tray	IBD-4030	
IV Pole	IBD-4040	
WARRANTY		
Standard Warranty	1 Year	Subject to MTTS policy

MTTS

MEDICAL TECHNOLOGY TRANSFER AND SERVICES Co., LTD
ISO 13485: 2016 Certified Company

For more information or to place an order please contact us at:

@ info@mtts-asia.com

www http://mtts-asia.com

+84 24 37 666 521

No 26, Lane 41 An Duong Vuong Street
Phu Thuong Ward, Tay Ho District
Hanoi, Vietnam

Distributor: