

Wallaby<sup>™</sup>  
warmer

LifeKit by **MTTS**



CE 2265

# WALLABY

An efficient, comfortable device to **prevent and treat hypothermia.**



## Overview

- REMOTE TEMPERATURE and object sensing providing full patient protection from overheating if standard skin sensor is disconnected or faulty
- Angle adjusted bed (up to 12 degrees)
- Adjustable overhead unit and tray for X-ray procedures

### 1. Temperature Control System

- Automatic (servo controlled)
- Manual
- Prewarm

### 2. Displays/Indicators

- bright LCD display to monitor patient temperature, treatment time, total usage time, heater power level and sensor connection
- built in APGAR score timer
- LCD display with diagnostic and troubleshooting information
- total system time counter (lifetime)

### 3. Safety

- Heater element failure
- Sensor failure
- Power failure
- Sensor off the patient

### 4. Optional Items

- IV Pole
- Storage Shelf
- Storage Drawer
- Cylinders holder



## Dimensions and Weight

|                             |                         |
|-----------------------------|-------------------------|
| Unit dimensions             | 110 cm x 78 cm x 174 cm |
| Shipping (boxed) dimensions | 38 cm x 122 cm x 82 cm  |
| Net Weight                  | 64 kg                   |
| Shipping Weight             | 78 kg                   |

## Performance Specifications

|                |                                                                     |
|----------------|---------------------------------------------------------------------|
| Application    | Treatment of neonatal and pediatric patients to prevent hypothermia |
| Heater Element | Type: Radiant<br>Power: 650W                                        |
| Procedure Lamp | Type: LED<br>Power: 3W                                              |

## Temperature

|                           |                          |
|---------------------------|--------------------------|
| Patient Temperature range | 34.0°C - 37.5°C in 0.1°C |
| Temperature measurement   | 10.0°C - 50.0°C in 0.1°C |
| Operating Temperature     | 18.0°C - 35.0°C          |

## Displays / Indicators

|                     |                    |
|---------------------|--------------------|
| Set temperature     | Total usage time   |
| Patient temperature | Heater power level |
| Temperature Alarm   | Sensor connected   |
| Treatment time      | Power indicator    |
| APGAR score timer   |                    |

## Extra Features

- Adjustable overhead unit for X-ray procedures
- Tray for X-ray films
- Adjustable bed angle (up to 12 degrees)
- Shelf for monitoring equipment

## Requirements

|                     |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Power Requirements  | 220-240VAC, 47/63Hz                                                                                              |
| Agency Requirements | ISO 13485: 2016<br>EN 60601-1:2006/AC2010<br>EN 60601-1-2:2007/AC2010<br>EN 60601-6:2010<br>IEC 60601-2-21: 2003 |

## Ordering Information and Servicing

| MAIN UNIT                         | ITEM ORDER NO |           |
|-----------------------------------|---------------|-----------|
| Warmer Unit                       | WLB-1000      | Main Unit |
| OPTIONAL ITEMS                    | ITEM ORDER NO |           |
| Skin Sensor - tape                | WLB-1011      |           |
| Sticky tape for skin sensor (100) | WLB-1013      |           |
| Mattress                          | WLB-1020      |           |
| IV Pole                           | WLB-1030      |           |
| Storage Shelf                     | WLB-1040      |           |
| Storage Drawer                    | WLB-1050      |           |
| Cylinders Holder                  | WLB-1060      |           |
| Baby Bonnet                       | WLB-1070      |           |
| Documents Holder                  | WLB-1080      |           |

## WARRANTY

|                   |        |                        |
|-------------------|--------|------------------------|
| Standard Warranty | 1 year | Subject to MTTs policy |
|-------------------|--------|------------------------|



MEDICAL TECHNOLOGY TRANSFER AND SERVICES Co., LTD  
ISO 13485: 2016 Certified Company

@ info@mtts-asia.com

www http://mtts-asia.com

+84 243 766 6521

No 26, Lane 41 An Duong Vuong Street  
Phu Thuong Ward, Tay Ho District  
Hanoi, Vietnam

Distributor: